Novel Coronavirus (SARS-COV-2) Antigen Rapid Test Cassette (Swab)

25 Tests/kit

Features:

Testing Time : 10 Mins

Sensitivity: 96.17%

Specificity: 99.9%


Product name	Ref.	Size	Storage Temp
Novel Coronavirus (SARS-COV-2) Antigen Rapid Test Cassette (Swab)	K511416D	25 Test/kit	2-30 °C
Novel Coronavirus (SARS-COV-2) Antigen Rapid Test Cassette (Swab)	K511416D	5 Test/kit	2-30 °C

REALy

CE

Realy Novel Coronavirus (SARS-COV-2) Antigen Rapid Test Cassette (Swab) is a rapid chromatographic immunoassay for the qualitative detection of specific antigens to SARS-CoV-2 present in human nasopharynx.

DIRECTIONS FOR USE:

Novel Coronavirus (SARS-CoV-2) Antigen Rapid Test Cassette (swab) Package Insert

A RAPID TEST FOR THE QUALITATIVE DETECTION OF NOVEL CORONAVIRUS ANTIGENS IN NASOPHARYNGEAL SWAB.

For professional In Vitro Diagnostic Use Only.

INTENDED USE

The Novel Coronavirus (SARS-CoV-2) Antigen Rapid Test Cassette (swab) is an *in vitro* diagnostic test for the qualitative detection of novel coronavirus antigens in Nasopharyngeal swab, using the rapid immunochromatographic method. The identification is based on the monoclonal antibodies specific for the novel coronavirus antigen. It will provide information for clinical doctors to prescribe correct medications.

SUMMARY

The Novel Coronavirus (SARS-CoV-2) Antigen Rapid Test Cassette (swab) is used for the *in vitro* qualitative detection of novel coronavirus in the throat swabs, sputum samples of suspected pneumonia patients infected by novel coronavirus, suspected clustering cases and others needing diagnosis or differential diagnosis for novel coronavirus.

The definitions of "suspected cases" and "patients with suspected aggregated cases" and other groups are implemented with reference to the "Diagnosis and Treatment Plan for Pneumonia Infected in novel coronavirus" and "Monitoring Plan for Pneumonia Infected in novel coronavirus" and other documents (current version) issued by CDC.

The product is only used for auxiliary diagnosis of related cases and emergency reserve for *in vitro* diagnosis of this epidemic during the pneumonia epidemic infected by novel coronavirus (SRAS-CoV-2) since December 2019 and it cannot be used as routine *in vitro* diagnostic reagents in clinical practice. The kit shall comply with the relevant requirements of the "Diagnosis and Treatment Plan for Pneumonia Infected in novel coronavirus" and "Prevention and Control Plan for Pneumonia Infected in novel coronavirus" and other documents in use.

The detection results of this kit are for clinical reference only and should not be used as the sole criteria for clinical diagnosis. It is recommended to conduct a comprehensive analysis on the condition in combination with the clinical manifestations and other laboratory tests.

PRINCIPLE

The Novel Coronavirus (SARS-CoV-2) Antigen Rapid Test Cassette (swab) is an immunochromatographic membrane assay that uses highly sensitive monoclonal antibodies to Novel coronavirus.

The test device is composed of the following three parts, namely sample pad, reagent pad and reaction membrane. The whole strip is fixed inside a plastic device. The reagent membrane contains the colloidal-gold conjugated with the monoclonal antibodies against Novel coronavirus; the reaction membrane contains the secondary antibodies for Novel coronavirus, and the polyclonal antibodies against the mouse globulin, which are pre-immobilized on the membrane.

When the sample is added into the sample window, conjugates dried in the reagent pad are dissolved and migrate along with the sample. If Novel coronavirus is present in the sample, a complex formed between the anti- Novel coronavirus conjugate and the virus will be caught by the specific anti- Novel coronavirus monoclonal coated on the T region.

Whether the sample contains the virus or not, the solution continues to migrate to encounter another reagent (an anti-mouse IgG antibody) that binds the remaining conjugates, thereby producing a red line on the region C.

REAGENTS

The reagent membrane contains the colloidal-gold conjugated with the monoclonal antibodies against Novel coronavirus; the reaction membrane contains the secondary antibodies for Novel coronavirus, and the polyclonal antibodies against the mouse globulin, which are pre-immobilized on the membrane.

PRECAUTIONS

- For *in vitro* diagnostic use only.
- Do not use after the expiration date.
- Ensure foil pouch containing test device is not damaged before opening for use.
- Perform test at room temperature 15 to 30°C.
- Wear gloves when handling the samples, avoid touching the reagent membrane and sample window.
- All samples and used accessories should be treated as infectious and discarded according to local regulations.
- Avoid using bloody samples.

STORAGE AND STABILITY

Store the Novel Coronavirus (SARS-CoV-2) Antigen Rapid Test Cassette (swab) at room temperature or refrigerated (2-30°C). Do not freeze. All reagents are stable until the expiration dates marked on their outer packaging and buffer vial.

SPECIMEN COLLECTION AND PREPARATION

1. Specimen collection:

It is applicable to the diagnosis of the Novel coronavirus from the samples of Nasopharyngeal swab. Use freshly collected samples for optimal test performance. Inadequate sample collection or improper sample handling may yield a false-negative result.

Nasopharyngeal Swab

For nasal swab completely insert the sterilized swab supplied in this kit into the nasal basin, and swab several times to collect the epidermal cells from the mucus.

For throat swab completely insert the sterilized swab supplied in this kit into the posterior pharynx, tonsils and other inflamed areas. Avoid touching the tongue, cheeks and teeth with the swab.

It is recommended to collect sample from Nasopharyngeal for more accurate results.

2. Specimen preparation:

1) Add 10 drops (about 0.3 ml) of the sample extraction buffer into the extraction tube supplied in this kit, and put it on the tube stand.

2) Nasopharyngeal Swabbing

Insert the swab into the extraction tube which contains 10 drops (about 0.3 ml) of the sample extraction buffer. Rotate the swab inside the tube using a circular motion to roll the side of the extraction tube so that liquid is expressed and reabsorbed from the swab. Remove the swab. The extracted solution will be used as test sample.

MATERIALS

Materials provided

- Test Device
- Sterilized Swab
- Extraction Tube
- Tube Stand
- Package Insert
- Nozzle With Filter
- Sample Extraction Buffer

Materials required but not provided

- Timer

DIRECTIONS FOR USE

Allow the test, specimen, extraction buffer to equilibrate to room temperature (15-30°C) prior to testing.

1.Remove the test device from the sealed foil pouch and use it as soon as possible. Place the test device on a clean and level surface. Best results will be obtained if the assay is performed immediately after opening the foil pouch.

2.Unscrew the whole cap of the specimen collection tube,


3.Add 10 drops (about 0.3 ml) of the sample extraction buffer into the extraction tube.

4.Place the sterilized swab specimen in the sample extraction buffer. Rotate the swab for approximately 10 seconds while pressing the head against the inside of the tube to release the antigen in the swab.

5.Remove the sterilized swab while squeezing the sterilized swab head against the inside of Buffer as you remove it to expel as much liquid as possible from the swab. Discard the sterilized swab in accordance with your biohazard waste disposal protocol.

6.Screw on and tighten the cap onto the specimen collection tube, then shake the specimen collection tube vigorously to mix the specimen and the sample extraction buffer.. See illustration 4.

7.Add 3 drops of the solution (approx.80ul) to the sample well and then start the timer. Read the result at 10-20 minutes. Do not interpret the result after 20 minutes.

INTERPRETATION OF RESULTS

(Please refer to the illustration above)

POSITIVE: Two red lines appear. One red line appears in the control region(C), and one red line in the test region(T). The shade of color may vary, but it should be considered positive whenever there is even a faint line.

NEGATIVE: Only one red line appears in the control region(C), and no line in the test region(T). The negative result indicates that there are no Novel coronavirus particles in the sample or the number of viral particles is below the detectable range.

INVALID: No red line appears in the control region(C). The test is invalid even if there is a line on test region(T). Insufficient sample volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the test procedure and repeat the test using a new test device. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

LIMITATIONS

- The Novel Coronavirus (SARS-CoV-2) Antigen Rapid Test Cassette (swab) is an acute-phase

screening test for qualitative detection. Sample collected may contain antigen titles below the reagent's sensitivity threshold, so a negative test result does not exclude infection with novel coronavirus.

• The Novel Coronavirus (SARS-CoV-2) Antigen Rapid Test Cassette (swab) detects viable and non-viable novel coronavirus antigen. Test performance depends on antigen load in the sample and may not correlate with cell culture performed on the same sample. A positive test does not rule out the possibility that other pathogens may be present. Therefore, the results must be compared with all other available clinical and laboratory information to make an accurate diagnosis.

• A negative test result may occur if the level of extracted antigen in a specimen is below the sensitivity of the test or if a poor quality specimen is obtained.

• Performance of the test has not been established for monitoring antiviral treatment of novel coronavirus.

• Positive test results do not rule out co-infections with other pathogens.

• Negative test results are not intended to rule in other coronavirus infection except the SARS-CoV-1

• Children tend to shed virus for longer periods of time than adults, which may result in differences in sensitivity between adults and children. List.

• A negative result may occur if the concentration of antigen or antibody in a specimen is below the detection limit of the test or if the specimen was collected or transported improperly, therefore a negative test result does not eliminate the possibility of SARS-CoV-2 infection, and should be confirmed by viral culture or an molecular assay or ELISA

PERFORMANCE CHARACTERISTICS

Clinical Evaluation

Clinical evaluation was performed to compare the results obtained by Novel Coronavirus (SARS-CoV-2) Antigen Rapid Test Cassette (swab) and PCR. The results were summarized below:

Table: Novel Coronavirus (SARS-CoV-2) Antigen Rapid Test Cassette (swab) vs. PCR

	SARS-CoV-2 Ag Rapid Test		Total Result
	+	-	
PCR	27	5	32
Total Results	27	200	232

Relative sensitivity: 84.4%

Relative specificity: >99%

Overall agreement: 97.8%

Cross Reaction

No cross reaction has been confirmed of the Novel Coronavirus (SARS-CoV-2) Antigen Rapid Test Cassette (swab) with the following pathogens:

①Bacteria

Acinetobacter baumannii, Bordetella pertussis, Branhamella catarrhalis, Candida albicans, Candida glabrata, Cardiobacterium hominis, Eikenella corrodens, Enterococcus faecalis, Enterococcus gallinarum, Escherichia coli, Group C streptococcus, Group G streptococcus, Haemophilus aphrophilus, Haemophilus influenzae, Haemophilus paraphilic, Klebsiella pneumoniae, Neisseria gonorrhoeae, Peptococcus asaccharolyticus, Peptostreptococcus anaerobic, Proteus mirabilis, Proteus vulgaris, Pseudomonas aeruginosa, Serratia marcescens, Staphylococcus epidermidis, Streptococcus agalactiae (group B), Streptococcus mutans, Streptococcus pneumoniae, Streptococcus pyogenes (group A), Veillonella parvula

②Virus

Influenza A, Influenza B, Adenovirus Type 1 ~ 8, 11, 19, 37, Coxsackie virus Type A16, B1 ~ 5, Cytomegalovirus, Echovirus Type 3, 6, 9, 11, 14, 18, 30, Enterovirus Type 71, HSV-1, Mumps virus, Typhoid simple herpes virus, Parainfluenza virus Type 1 ~ 3, Poliovirus Type 1 ~ 3, Respiratory syncytial virus, Rhinovirus Type 1A, 13, 14, Type I simple herpes virus, MERS.

③Mycoplasma etc.

No cross reaction with Chlamydia pneumoniae, Chlamydia psittaci, Chlamydia trachomatis, Mycoplasma pneumoniae.


SYMBOLS

Symbol	Meaning	Symbol	Meaning
IVD	In vitro diagnostic medical device		Storage temperature limit
		EC REP	Authorized representative in the European Community
	Date of Manufacture		Use by date
	Do not reuse		Consult instruction for use
LOT	Batch code	CE	Meet the requirements of EC Directive 98/79/EC

HANGZHOU REALY TECH CO., LTD.

4th Floor, #12 Building, Eastern Medicine Town, Xiasha Economic & Technology Development, 310018 Hangzhou, Zhejiang, P. R. China

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Manufacturer:

Name: HANGZHOU REALY TECH CO., LTD.

Address: 4th Floor, #12 Building, Eastern Medicine Town, Xiasha Economic & Technology Development, 310018 Hangzhou, Zhejiang, P.R.China

Product/s	Catalogue number
Novel Coronavirus (SARS-CoV-2) Antigen Rapid Test Cassette (swab)	K511416D

Category: Other Devices (All devices except Annex II and self-testing devices)

Conformity assessment route: Annex III, except Point 6, of Directive

Applicable Standards: EN ISO 13485:2016; EN ISO 15223-1:2016;
EN ISO 14971:2012; EN ISO 13612:2002; EN ISO 17511:2003;
EN ISO 18113-1:2011; EN ISO 18113-2:2011; EN ISO 23640:2015.

We, the Manufacturer, herewith declare with sole responsibility that our product/s mentioned above meet/s the provisions of the Directive 98/79/EC of the European Parliament and of the Council on In-Vitro Diagnostic Medical Devices.

We hereby explicitly appoint Luxus Lebenswelt GmbH, located at Kochstr. 1, 47877, Willich, Germany to act as our European Authorised Representative as defined in the aforementioned Directive.

Hangzhou 2020.8.17

(Place and date of issue)

(Signature and position)

Signed for and on behalf of the manufacturer

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Formblatt für In-vitro-Diagnostika / Form for In Vitro Diagnostic Medical Devices

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Sample Name: Novel Coronavirus (SARS-CoV-2) Antigen Rapid Test Cassette (swab)

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Shanghai Research Institute of Chemical Industry Testing Co., Ltd

海 运
By Sea

海 运
By Sea

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Certification for Safe Transport of Chemical Goods


NO. 2020088765

Page 1 / 2

样品名称 Sample Name	中文 Chinese	新型冠状病毒抗原检测试剂盒
	英文 English	Novel Coronavirus (SARS-CoV-2) Antigen Rapid Test Cassette (swab)
委托单位 Consignor	杭州睿丽科技有限公司	
生产单位 Manufacturer	杭州睿丽科技有限公司	
检验方法、程序 Inspection Methods and Procedures	国际海事组织《国际海运危险货物规则》(2018版) IMO International Maritime Dangerous Goods Code (2018 Edition)	
样品外观与气味 Appearance & Odor	多种颜色纸盒(内含白色塑料测试板及无色透明液体), 稍有气味 Multicolor Paper box(containing white plastic test board and colorless transparent liquid), Weak odor	
IDENTIFICATION CONCLUSION	1. 危险性识别(Hazards identification) 无。 None. 2. 海运按照IMO IMDG Code办理的类项(Suggestion according to IMO IMDG Code) 根据特殊规定286, 可按非限制性货物条件办理。 The substance is not subject to IMO IMDG Code according to special provision 286. 3. 包装要求(Packaging requirements) 无。 None.	
备注 Comment	无。 None.	

检验日期: 2020-08-31 签发日期: 2020-08-31 生效日期: 2020-08-31
Inspection Date: 2020-08-31 Issue Date: 2020-08-31 Effective Date: 2020-08-31

批准
Approver: 

审核
Checker: 

主检
Appraiser: 

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Certification for Safe Transport of Chemical Goods

NO. 2020088765

Page 2 / 2

鉴定项目 Identification Items	鉴定结果 Identification Conclusion Results
爆炸危险性鉴定 Identification of Explosive Hazard	该货物不属于爆炸品。 The product is not classified in Explosives.
易燃危险性鉴定 Identification of Flammable Hazards	该货物不属易燃危险品。 The product is not classified in flammable substance.
氧化危险性鉴定 Identification of Oxidative Hazards	该货物不属于氧化剂和有机过氧化物。 The product is not classified in oxidizing substances and organic peroxides.
毒害及传染危险性鉴定 Identification of Toxic & Infectious Hazards	该货物不属于有毒和感染性物质。 The product is not classified in toxic and infectious substances.
放射危险性鉴定 Identification of Radioactive Hazard	该货物无放射危险性。 The product is not classified in radioactive material.
腐蚀危险性鉴定 Identification of Corrosive Hazard	该货物不属于腐蚀品。 The product is not classified in corrosives.
其他危险性鉴定 Identification of other Hazards	该货物无其它危险性。 The product presents no other dangerous properties.

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空运
By Air

NO.2020088764

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Certification
for Safe Transport of Chemical Goods

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Shanghai Research Institute of Chemical Industry Testing Co., Ltd

空运
By Air

空运
By Air

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Certification for Safe Transport of Chemical Goods

NO. 2020088764

Page 1 / 2

样品名称 Sample Name	中文 Chinese	新型冠状病毒抗原检测试剂盒
	英文 English	Novel Coronavirus (SARS-CoV-2) Antigen Rapid Test Cassette (swab)
委托单位 Consignor	杭州睿丽科技有限公司	
生产单位 Manufacturer	杭州睿丽科技有限公司	
检验方法、程序 Inspection Methods and Procedures	国际航空运输协会《危险品规则》61版 IATA Dangerous Goods Regulations (DGR) 61st Edition	
样品外观与气味 Appearance & Odor	多种颜色纸盒(内含白色塑料测试板及无色透明液体), 稍有气味 Multicolor Paper box(containing white plastic test board and colorless transparent liquid), Weak odor	
IDENTIFICATION 鉴定结论 CONCLUSION	<p>1. 危险性识别 (Hazards identification) 无。 None.</p> <p>2. 空运按照IATA DGR办理的类项 (Suggestion according to IATA DGR) 根据特殊规定A122, 可按非限制性货物条件办理。 The substance is not subject to IATA DGR according to special provision A122.</p> <p>3. 包装要求 (Packaging requirements) 无。 None.</p>	
备注 Comment	无。 None.	

检验日期: 2020-08-31 签发日期: 2020-08-31 生效日期: 2020-08-31
Inspection Date: 2020-08-31 Issue Date: 2020-08-31 Effective Date: 2020-08-31

批准
Approver: 

审核
Checker: 

主检
Appraiser: 

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Certification for Safe Transport of Chemical Goods

NO. 2020088764

Page 2 / 2

鉴定项目 Identification Items	鉴定结果 Identification Conclusion Results
爆炸危险性鉴定 Identification of Explosive Hazard	该货物不属于爆炸品。 The product is not classified in Explosives.
易燃危险性鉴定 Identification of Flammable Hazards	该货物不属易燃危险品。 The product is not classified in flammable substance.
氧化危险性鉴定 Identification of Oxidative Hazards	该货物不属于氧化剂和有机过氧化物。 The product is not classified in oxidizing substances and organic peroxides.
毒害及传染危险性鉴定 Identification of Toxic & Infectious Hazards	该货物不属于有毒和感染性物质。 The product is not classified in toxic and infectious substances.
放射危险性鉴定 Identification of Radioactive Hazard	该货物无放射危险性。 The product is not classified in radioactive material.
腐蚀危险性鉴定 Identification of Corrosive Hazard	该货物不属于腐蚀品。 The product is not classified in corrosives.
其他危险性鉴定 Identification of other Hazards	该货物无其它危险性。 The product presents no other dangerous properties.

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25 tests/kit

500 tests/carton

Carton size: 45*44*28cm

Volume: 0.056CBM

Gross weight per carton: 7.5KG

Volume weight per carton via air cargo: 9.5KG

Volume weight per carton via express: 11.5KG

5 tests/kit

1075 tests/carton

Carton size: 65*55*45cm

Volume: 0.16CBM

Gross weight per carton: 20KG

Volume weight per carton via air cargo: 27KG

Volume weight per carton via express: 32KG